FLUORESCENT MOLECULAR ROTORS FOR MEASURING INTRACELLULAR VISCOSITY

Anna Blood¹, Nopondo Esemoto², Adam Meares², Miji Jeon², Marcin Ptaszek², Songon An² ¹Department of Chemistry and Biochemistry, New College of Florida, 5800 Bay Shore Rd, Sarasota, FL 34243 ²Department of Chemistry and Biochemistry, University of Maryland, Baltimore County, 1000 Hilltop Circle, Baltimore, MD 21250

Intro	Synthesis	Result
A method to visualize intracellular viscosity will allow us to develop our understanding of assembly mechanisms of biological pathways in single cells. When proteins interact, colocalization of enzymes form clusters. ¹	The PEG chains were synthesized following the procedure from the Li group (Scheme 1). ⁴ The core of target compound BDP 2 was synthesized following the procedure of preexisting compound BDP 1 (Scheme 2). ²	Human breast carcinoma Hs578T cells were treated with BDP 1 and BDP 2 and imaged by DIC microscopy and fluorescence microscopy.
which results in a more viscous environment. Visualization can inform us on which proteins are interacting and where the interactions are taking place. The goal of my research is to develop a novel type of fluorescent	The PEG chains were connected to the core by a click reaction (Scheme 2).	BDP 1 forms highly fluorescent clusters in the cell, which appear directly next to vesicles imaged with DIC microscopy. BDP 2 also forms clusters in the cell but the fluorescence intensity is much weaker. Although this is

environments, rotation quenches fluorescence, whereas in more viscous environments, limited rotation results in more intense fluorescence.

Background



Less Viscous More Rotation Less Florescent

More Viscous Less Rotation More Florescent

When a molecule absorbs energy and moves to an excited energy state, it can relax back down to the ground state by releasing photons, which results in fluorescence, or by non-radiative internal conversions by rotating around the two conjugated portions of the molecule.

In less viscous media, the molecule has more freedom to rotate, so more energy will be lost by non-radiative internal conversion. Less energy is available to be released as photons in order for the molecule to relax back down to the ground state. So, the rotor will be less fluorescent in less viscous media.



Reaction Conditions: (1a) TsCl, NaOH, THF/H₂O (1:1), RT, 3 h, 65%; (1b) NaN₃, DMF, 90°C, overnight, 37%; (2a) DDQ, TEA, BF₃•OEt₂, RT, 1 h, 35%; (2b) AlCl₃, THF, Propargyl Alcohol, 40°C, 3 h, RT, 1 h, 78%; (2c) Tetra PEG Azide, CuSO₄•5H₂O, Ascorbic Acid, Acetone Water (5:1), Microwave, 150 W, 65°C, 15 m, 24%.

Spectroscopic Characterization

the dye in the cell, and the longer fluorescence exposure time results in increased photobleaching. Additionally, the pattern of localization of BDP 2 is not as clear as that of BDP 1, but it is hypothesized that BDP 2 is accumulated in the Golgi apparatus.



Figure 3. Images of BDP 1 [A-C], BDP 2 [D-F], and BDP 2 fixed [G-I] in Hs578T cells. DIC microscopy [left], fluorescence microscopy [middle], merged [DIC (red), fluorescence (green), right]

In a more viscous environment, rotation will be inhibited, so less energy will be lost to non-radiative internal conversion. More photons will be released in order for the molecule to relax back down to the ground state. So, in more viscous media, the rotor will have stronger fluorescence.

Molecular Design



The molecular rotor is based upon boron-dipyrromethene (BODIPY), a commonly used green fluorescent dye. BDP 1 was used as a model for target compound BDP 2.² While the BODIPY portion of the molecule is

The spectroscopic characterization of the molecular rotor shows that the quantum yields increase as solvent viscosity increases, as anticipated. Additionally, the rotor is soluble and fluorescent in water, which is important for use in biological applications, with an absorption maximum at 507 nm and an emission maximum at 531 nm.

Solvent	Dynamic Viscosity at 25°C [cP]	Fluorescence Quantum Yield [Φ _F]
Toluene	0.56	0.03
DI Water	0.89	0.15
PEG 400	90.0	0.21

Figure 1. Summary of BDP 2 quantum yield measurements in toluene, DI water, and PEG 400.



Conclusion and Future Direction

In summary, a novel fluorescent molecular rotor was synthesized to visualize intracellular viscosity. In less viscous environments, rotation quenches fluorescence, whereas in more viscous environments, limited rotation results in stronger fluorescence. This was supported by quantum yield measurements in vitro, which increased as solvent viscosity increased. Imaging in Hs578T cells indicated that the compound requires optimization for better visualization in cells.

The compound will be optimized by changing the substituent to observe how fluorescence intensity and intracellular localization are affected. By transfecting the cells with glucosome proteins, it can be determined if the localization of the fluorescent dye is related to the enzyme clusters and provide insight on the assembly mechanism of the glucosome proteins.⁵ Additionally, BDP 2 will be imaged with the Golgi apparatus to be evaluated for its localization in living cells.



hydrophobic, the polyethylene glycol (PEG) chains are hydrophilic, making the compound amphiphlic for localization in the cytoplasm. The PEG chains were chosen due to their lack of charge and known biocompatibility.³ The BODIPY meso aryl substituent was chosen to enable free rotation in order to quench fluorescence. BDP 1 has methyl groups on the pyrrole rings, which inhibits the rotation of the phenyl group. BDP 2 was synthesized without the methyl groups so rotation would not be inhibited.

Figure 2. Summary of BDP 2 absorption and emission spectra in DI water.

We thank the National Science Foundation Research Experiences for Undergraduates (NSF-REU) Sites program. This material is based upon work supported by the National Science Foundation under CHE-1460653.